

# Development of a Total Beryllium Field Analyzer

Amy Ekechukwu, Paula Cable-Dunlap, Scott McWhorter, Charles Shick, Jr., Jeff Degange,

# Executive Summary

- **SRS assembled cross-division team to address Be analyses and digestion**
- **Need for Be field analysis compliant with 10 CFR 850 (real-time and near-real-time)**
- **SRS Field deployable analyzer - combines digestion method with colorimetric analysis**

# Why is Beryllium Important?

- **Beryllium inhalation can result in development of Chronic Beryllium Disease (CBD)**
- **Several DOE sites have beryllium melting, casting, grinding, and machining operations (past or current)**
- **The DOE Beryllium Rule (10 CFR Part 850) established free-release limit of removable Be surface contamination at 0.2 micrograms per 100 cm<sup>2</sup>**

# Be Analysis: Standard Protocol

- **Sample Digestion**
  - Nitric acid (7300); nitric/sulfuric (7102)
- **Analysis**
  - Usually Inductively coupled plasma emission spectrometric (ICP-ES) analysis (7300)
  - Mass Spec (ICP-MS) , Atomic Absorbtion (GFAA) (7102) also used
- **Laboratory analyses - turnaround time days to weeks**
- **Need for real-time or near-real-time method**

# Be Field Analyzer Development

- **SRS Be Team addressing method improvements**
  - robust sample digestion for BeO
  - alternative analytical methods
    - colorimetric technique
- **Field Analyzer - combines digestion method and colorimetric method into field portable instrument**

# Requirements for Field Analyzer

- **Able to measure all forms of Be (oxide, metal)**
- **Meet free-release detection limit (0.2  $\mu\text{g}/\text{swipe}$ )**
- **Portable**
- **Results at or near real-time**
- **Easy to manage wastes**
- **Simple and rugged**

# SRS Field Analyzer

- **Combines digestion and colorimetric method**
- **Cost-effective**
  - 75% cost reduction estimate
- **Relatively short turnaround time**
  - hours versus days
- **Provides low detection limits required in DOE Be Rule**
- **Measures all forms of beryllium**

# Sample Digestion - Why?

- **Formation of highly insoluble BeO can occur in DOE Be processes**
- **Standard NIOSH method (7300) does not address oxide (BeO)**
- **For “total” beryllium measurement, all Be must be in solution**
- **Sample matrix - Ghostwipe™ swipes (Environmental Express)**



# Typical Ghostwipe Swipe



# ICP-MS Analysis of Machine Shop Swipe Matrix

## Swipe Matrix Composition

Element	SQ Standard Known Concn. (ppb)	SQ Standard Measured Concn. (ppb)	Blank Swipe (µg/swipe)	Swipe Matrix (µg/swipe)
Li		0.0038	0.016	0.93
Be	100	100	<0.0004	0.014
B		0.24	1.2	34
Na		3	600	940
Mg	100	100	125	310
Al		1.9	3.5	16000
Sc		0.11	0.20	1.9
Ti		0.16	0.33	21
V		0.018	0.032	120
Cr		0.25	0.6	10000
Mn		0.24	0.29	1000
Fe		41	11	51000
Co	100	97	0.024	250
Ni	100	100	0.55	3400
Cu		0.9	1.8	7500
Zn		0.28	9.5	170
Ga		0.017	0.035	6.6
Sr		0.013	0.14	7.3
Y		0.005	0.0018	0.21
Zr		0.0098	0.05	1.4
Nb		<0.002	0.0018	7.4
Mo		0.1	0.15	460
Ru		<0.0098	<0.0005	4
Pd		<0.0097	0.0043	4.2
Ag		<0.004	0.027	3.5
Cd		<0.022	0.26	2.1
In	100	100	0.39	13
Sn		0.02	0.13	45
Sb		<0.0081	0.07	3.5
Ba		0.025	0.46	29
Ce	100	100	0.0047	0.76
Ta		<0.001	0.0017	0.63
W		0.45	0.39	520
Re		<0.0028	0.015	0.43
Au		<0.002	0.0011	0.13
Pb	100	100	1.3	55
Bi	100	100	0.018	1.5
U	100	100	0.016	0.14


## Key Matrix Components

Element	SQ Standard Known Concn. (ppb)	SQ Standard Measured Concn. (ppb)	Blank Swipe (µg/swipe)	Swipe Matrix (µg/swipe)
Be	100	100	<0.0004	0.014
Mg	100	100	125	310
Al		1.9	3.5	16000
V		0.018	0.032	120
Cr		0.25	0.6	10000
Mn		0.24	0.29	1000
Fe		41	11	51000
Co	100	97	0.024	250
Ni	100	100	0.55	3400
Cu		0.9	1.8	7500
Zn		0.28	9.5	170
Mo		0.1	0.15	460
In	100	100	0.39	13
Ce	100	100	0.0047	0.76
W		0.45	0.39	520
Pb	100	100	1.3	55
U	100	100	0.016	0.14

# Digestion Methods Evaluated

- Multiple digestion protocols tested
- Fired Be-oxide particularly difficult to digest
- Most promising digestion involves  $\text{HNO}_3$ ,  $\text{HCl}$ ,  $\text{HF}$ ,  $\text{H}_2\text{O}_2$ , Heating

# EPA Modified 3050B Digestion



# Colorimetric Method - 3 Stages

- **SRTC modified LANL colorimetric wipe method**
- **Developed semi-quantitative liquid technique**
- **Used spectroscopy, enhanced chemistry to obtain linear results at required detection limit**

# Stage 1 - Modified Swipe Method

- **LANL Method - Treated swipes with reagents, wiped surface, observe color**
- **SRS - Used same chemistry**
  - wiped surface with clean wipe
  - add EDTA - Complex metals
  - buffer to pH 9-10 - Make Be available for reagent
  - add color reagent, Chromium Azurol S (CAS) and observed color change
- **Could not meet free-release criteria**

## Stage 2 - Colorimetric Method


- **Swipe surface**
- **Wash swipe with acid to solublize Be**
- **Perform chemistry on acid wash of swipe**
- **Observe color of liquid - positive visible result at 1 ug/swipe**
- **CAS (red) interfered with visible detection of CAS-Be (purple) complex at necessary detection limit (0.2 ug/swipe)**

# Stage 3 - Absorbance Method

- **Needed to lower detection limits to 0.2 ug/swipe**
- **Used spectrophotometer**
  - wash Be off of swipe and react with colorimetric reagent
  - measure absorbance.
- **Addition of complexing agent to shift CAS-Be absorbance away from CAS absorbance**




# Initial Calibration Plot




11-Feb-04

# Be-CAS-CTAB Optimized Calibration




# On-Line Swipe Digestion and Colorimetry (OSDAC) Interface

## Proof of Concept Device



- On-line swipe digestion
- Automated sample transfer to dilution/mixing station then to analyzer
- Separate control for acid flush of transfer lines to avoid sample cross-contamination

## Fieldable OSDAC



- Fully automated, integrated swipe digestion, sample transfer and analysis
- User friendly operation
- Robust, compact, lightweight
- Modular to prevent sample cross contamination

# Planned Tasks

- **Test digestion method on high fired Be-oxide - Fired Be-oxide being prepared for testing**
- **Decrease digestion time**
- **Test colorimetry on matrix matched samples containing interfering metals**
- **Test field analyzer using digestion and colorimetric chemistry**

# Path Forward

Develop Field  
BeO Digestion  
Method (in progress)

Develop Field  
Colorimetric Method  
(in progress)

Interface to  
Produce Field  
Analyzer  
(fabricated,  
awaiting  
chemistry)

Field Test (future)

Automate (future)

*In-situ* air monitor  
for Be (under development)

NIOSH  
Acceptance (future)

# Digestion Methods Under Evaluation

- Y-12 -  $\text{H}_2\text{O}_2$ , Conc.  $\text{H}_2\text{SO}_4$ , microwave
- 3050B -  $\text{HNO}_3$ ,  $\text{H}_2\text{O}_2$ ,  $\text{H}_2\text{SO}_4$
- 3050B Modified -  $\text{HNO}_3$ ,  $\text{HCl}$ ,  $\text{HF}$ ,  $\text{H}_2\text{O}_2$ , Dilute  $\text{HCl}$
- $\text{HNO}_3$ , sonication
- $\text{HNO}_3$ ,  $\text{H}_2\text{SO}_4$ ,  $\text{HF}$
- $\text{HCl}$ ,  $\text{CuSO}_4$ , sonication, heat (10 min.), 3:1  $\text{HNO}_3$ : $\text{H}_2\text{SO}_4$ , 2:1  $\text{HNO}_3$ :  $\text{HF}$

# Acknowledgements

Mike Brisson (SRS)  
Ed Kahal (SRS)  
Burney Hook (SRS)  
Scott McWhorter (SRS)  
Lori Chandler (SRS)  
Tony Woltermann (SRS)  
Alex Guanlao (SRS)  
June Hart (SRS)  
Art Jurgensen (SRS)  
Kathy Creek (LANL)

Chuck Coleman (SRS)  
Frank Pennebaker (SRS)  
Damon Click (SRS)  
Tom Oatts (Y-12)  
Steven Jahn (SRS)  
Steve LaMont (SRS)  
Tammy Taylor (LANL)  
Mike McCauley (McCauley  
Consultants)  
Ross Breese – AWE  
Aldermaston, U. K